


EcoFlow



NDC Accessories

| | |
|------------|----------------------|
| TM-NDCT013 | Screw On Cannula Cap |
|------------|----------------------|

Disposable Screw on Cap (sold separately)

- Packaged sterile
- Easy twist on and off
- Removable for easy insertion of grafts
- Universal for all sizes

NDC Cannula

| | |
|---------------|------------------------------|
| TM-NDCT557-MX | 5.5 mm x 72 mm Metal Cannula |
| TM-NDCT707-MX | 7.0 mm x 72 mm Metal Cannula |
| TM-NDCT857-MX | 8.5 mm x 72 mm Metal Cannula |

- Reduces waste in the operating room
- Windows at distal tip for visualization
- Unique thread profile for better retention
- Robust, metal obturator - Available non-cannulated or cannulated for going over a switching stick

NDC Obturator

| | |
|--------------|--|
| TM-NDCT055-C | 5.5 mm x 72 mm Cannula Obturator - Cannulated w/ Handle |
| TM-NDCT055-X | 5.5 mm x 72 mm Cannula Obturator - NonCannulated w/ Handle |
| TM-NDCT070-C | 7.0 mm x 72 mm Cannula Obturator - Cannulated w/ Handle |
| TM-NDCT070-X | 7.0 mm x 72 mm Cannula Obturator - NonCannulated w/ Handle |
| TM-NDCT085-C | 8.5 mm x 72 mm Cannula Obturator - Cannulated w/ Handle |
| TM-NDCT085-X | 8.5 mm x 72 mm Cannula Obturator - NonCannulated w/ Handle |

