




TIGON MEDICAL

GreenRepair[®]

Suture and Tape Selections

TigoCable, TigoTape, TigoLoop, TigoPass




Ordering Catalog


office@tigonmedical.com

Comprehensive Suture and Tape Selections

Tigon Medical offers a comprehensive selection of sutures and tapes with a variety of configurations including specialty options. The options are available in multiple color patterns, various width and lengths, and offered both with and without needles.



Softer and Less Abrasive





Testing was performed to determine how much force it would take for suture to pull through a tendon model. The competitor (Arthrex) tore through the model at ~30 N, where TigoCable pulled through at ~40 N and TigoTape at ~45 N.



**Save your fingers with
Tigon Medical's suture selection!**




Competitive Strength



Low Profile Knot




Needle Options



TigoLoop


TigoLoops are a needled, UHMWPE non-absorbable continuous suture loop that can be used in a variety of applications. Available in multiple configurations, tape or cable, and in various colors.



Part Number	Description
TST-154NL	1.5 mm White/Blue 26 mm Needle
TST-156NL	1.5 mm White/Black 26 mm Needle
TST-250NL	2.0 mm White 45 mm Needle
TST-258NL	2.0 mm White/Blue 45 mm Needle
TST-261NL	2.0 mm White/Black 45 mm Needle
TST-269NL	2.0 mm Blue/Black 45 mm Needle
TSS-504CNL	USP 5 White 45 mm Needle
TSS-505CNL	USP 5 White/Blue 45 mm Needle
TSS-509CNL	USP 5 White/Blue/Black 45 mm Needle

TigoPass


The TigoPass can be used to shuttle multiple sutures at once using the Dragon suture passer and allow for a luggage tag or cinch stitch technique. Currently available in cable and can be made with various colors.



Part Number	Description
TSS-209CL	USP 2 White/Black Non-Needled

TigoCable

Tigon Suture is a non-absorbable sterile suture made of ultra high molecular weight polyethylene (UHMWPE).




Available with or without suture needle and in various colors.

Size

- USP 2-0
- USP 2
- USP 5

Part Number	Description	
TSS-029CN	USP 2-0 White/Black 26 mm Needle	A
TSS-200C	USP 2 White Non-Needled	B
TSS-206C	USP 2 White/Blue No Needle	C
TSS-209CN	USP 2 White/Black 40 mm Needle	D
TSS-508CN	USP 5 White/Black 45 mm Needle	E



Part Number	Description	
1 TST-154	1.5 mm White/Blue Non-Needled	
2 TST-156	1.5 mm White/Black Non-Needled	
TST-168N	1.5 mm White/Blue 26 mm Needle	
TST-169N	1.5 mm White/Black 26 mm Needle	
3 TST-252N	2.0 mm White 45 mm Needle	
4 TST-254N	2.0 mm White/Blue 45 mm Needle	
5 TST-256N	2.0 mm White/Black 45 mm Needle	

TigoTape

Tigon Suture Tapes are an ultra high molecular weight polyethylene (UHMWPE) non-absorbable suture tape. The suture has round ends and a flat middle, improving soft tissue compression and force distribution. Available with or without suture needle and in various colors.

Size

- 1.5 mm
- 2.0 mm